

# Measuring Dihydrotestosterone (DHT) in Blood Serum for Research Purposes using Derivatization and LC-MS/MS

Debadeep Bhattacharyya<sup>1</sup>, Brad Hart<sup>1</sup>, Joe Di Bussolo<sup>1</sup>, Raidiri Castillo<sup>2</sup>, Ali Mustafa<sup>2</sup>, Catherine Lintag<sup>2</sup> and Hashim Othman<sup>2</sup>; <sup>1</sup>Thermo Fisher Scientific, San Jose, CA, <sup>2</sup>Bio-Reference Laboratories, Elmwood Park, NJ.

## ABSTRACT

We developed an LC-MS/MS method for researchers to accurately measure dihydrotestosterone (DHT) in blood serum from 25 to 2,000 pg/mL with a throughput of at least 12 injections per hour. This was done by liquid-liquid extraction of blood serum followed by derivatization with hydroxylamine and quantitative analysis using a 4-channel high-performance liquid chromatography (LC) system coupled to a tandem mass spectrometer (MS/MS). Oxime derivatives of DHT and its D<sub>3</sub> internal standard were separated from sample matrix components by gradient elution through an HPLC column containing solid-core silica particles with a polar-embedded reversed-phase bonded on its surfaces. Analytes were eluted to a heated electro-spray ionization probe of a triple-quadrupole mass spectrometer used for selected-reaction monitoring of the analytes. The desired analytical range of 25 to 2,000 pg/mL was achieved with inter- and intra-batch reproducibility of less than 7%, carryover less than 0.2% and acceptable correlation of specimen results with those from a reference laboratory. Throughputs of 12 to 24 injections/hour were achieved using one or two LC channels, respectively.

## INTRODUCTION

Dihydrotestosterone (DHT), also known as 5 $\alpha$ -dihydrotestosterone, androstanolone or stanolone, is an endogenous steroid hormone and a more potent androgen than testosterone. DHT is made by reduction of testosterone by 5 $\alpha$ -reductase. Researchers studying the physiology of DHT and control of its biosynthesis using 5 $\alpha$ -reductase inhibitors need to quantify this steroid within a range of 25 to 2,000 pg/mL (86.2 to 6,900 pmol/L) in blood serum. Since DHT does not ionize well by either atmospheric-pressure chemical ionization (APCI) or electrospray ionization (ESI), derivatization with hydroxylamine prior to liquid chromatography coupled to tandem mass spectrometry (LC-MS/MS) was necessary to reliably achieve the desired measuring range (1). Chemical structures and formulas of DHT and DHT-D3 internal standard are shown in Figure 1 along with mass-to-charge ratios (m/z) of their oxime derivatives. Fragmentation of DHT oxime is shown in Figure 2.

Figure 1. Analyte & internal standard (IS)





Figure 2. Fragmentation of DHT oxime



## MATERIALS AND METHODS


### Sample preparation

200  $\mu$ L aliquots of donor blood serum specimens, calibrators and quality control specimens (QCs), were spiked with DHT-D3 internal standard (IS) before being subjected to liquid-liquid extraction (LLE) with 2 mL methyl t-butyl ether (MTBE). After drying the MTBE extracts with heated nitrogen flow, the residue of each was reacted with hydroxylamine to form positive-ion oxime derivatives. After drying the derivative reagent, the residue of each was reconstituted with water and methanol (1:1) to a total volume of 200  $\mu$ L. Then, 50  $\mu$ L injections of each were made into the multi-channel LC-MS/MS system.

### LC-MS/MS methods


Using one or more channels of a Thermo Scientific™ Transcend™ LX-4 HPLC system, chromatographic separation of the steroid oximes from unwanted sample components was achieved by gradient elution through a Thermo Scientific™ Accucore™ aQ column (2.6  $\mu$ m, 100 x 2.1 mm), which was heated to 40°C. Figure 3 describes chromatographic conditions.

Figure 3. Chromatographic conditions for DHT-oxime



The Thermo Scientific™ TSQ Endura™ triple-quadrupole mass spectrometer was used with a heated electro-spray ionization (HESI) probe to selectively monitor analytes. Selected-reaction monitoring (SRM) of two transitions (quantitation and conformation) for DHT and DHT-D3 occurred within the 2-minute data window of the 4.5-minute LC method. Ion source and MS/MS conditions are described in Figure 4. The "Electro-clean" step is described in reference 2 and helped minimize charging effects.

Figure 4. Tandem mass spectrometry conditions



### Instrument control & data analysis


Thermo Scientific™ TraceFinder™ configured with Aria™ MX and TSQ Endura™ control software was used to control the LC-MS/MS systems, submit batches to the desired channels as well as to analyze data and report results.

## RESULTS

### Quantitation reliability

Typical DHT results for calibrators, QCs and specimens are shown in Figure 5. The desired measuring range from 25 to 2,000 pg/mL was achieved and was consistently linear ( $R^2 \geq 0.995$  with 1/X weighting). Carryover was less than 0.2%. Ion ratios (confirming/quant peak areas) averaged 54% for DHT-oxime and 46% for DHT-D3-oxime (IS). IS peak areas among calibrators & QCs averaged 19,200 with RSDs of 12.5%. IS peak areas from donor serum specimens ranged from 6,610 to 15,500 showing an average recovery of 59%, which indicated moderate ion-suppression by matrix. However, the IS in each sample preparation adequately compensated for such matrix effects.

Figure 5. Typical DHT quantitation results



As shown in Tables 1a and 1b, intra- and inter-batch precisions among 20 replicate injections from three pools (low, medium and high DHT levels) were less than 7% and 6% coefficient of variation (CV), respectively

Table 1. Intra- and inter-batch precisions

| a. Intra-batch results (pg/mL) |       |        |        |  | b. Inter-batch results (pg/mL) |       |        |        |      |
|--------------------------------|-------|--------|--------|--|--------------------------------|-------|--------|--------|------|
| Run                            | Low   | Medium | High   |  | Run                            | Low   | Medium | High   |      |
| 1                              | 172   | 590    | 1250   |  | Day 1                          | 1     | 185    | 520    | 1140 |
| 2                              | 182   | 550    | 1280   |  | 2                              | 208   | 500    | 1150   |      |
| 3                              | 193   | 580    | 1390   |  | 3                              | 197   | 560    | 1150   |      |
| 4                              | 190   | 590    | 1300   |  | 4                              | 195   | 580    | 1250   |      |
| 5                              | 179   | 570    | 1200   |  | 5                              | 204   | 540    | 1220   |      |
| 6                              | 177   | 570    | 1240   |  | Day 2                          | 1     | 195    | 550    | 1150 |
| 7                              | 187   | 610    | 1210   |  | 2                              | 207   | 550    | 1150   |      |
| 8                              | 174   | 600    | 1230   |  | 3                              | 182   | 570    | 1190   |      |
| 9                              | 188   | 560    | 1250   |  | 4                              | 199   | 550    | 1170   |      |
| 10                             | 171   | 570    | 1190   |  | 5                              | 188   | 580    | 1150   |      |
| 11                             | 185   | 570    | 1300   |  | Day 3                          | 1     | 201    | 520    | 1140 |
| 12                             | 160   | 540    | 1190   |  | 2                              | 192   | 530    | 1100   |      |
| 13                             | 163   | 600    | 1190   |  | 3                              | 186   | 490    | 1190   |      |
| 14                             | 195   | 520    | 1320   |  | 4                              | 202   | 530    | 1070   |      |
| 15                             | 163   | 590    | 1260   |  | 5                              | 193   | 500    | 1130   |      |
| 16                             | 173   | 550    | 1300   |  | Day 4                          | 1     | 209    | 520    | 1070 |
| 17                             | 175   | 590    | 1320   |  | 2                              | 201   | 540    | 1070   |      |
| 18                             | 165   | 580    | 1180   |  | 3                              | 191   | 510    | 1040   |      |
| 19                             | 160   | 550    | 1130   |  | 4                              | 189   | 560    | 1040   |      |
| 20                             | 184   | 550    | 1300   |  | 5                              | 192   | 510    | 1040   |      |
| Mean:                          | 176.8 | 571.5  | 1260.5 |  | Mean:                          | 195.8 | 535.5  | 1130.5 |      |
| StdDev:                        | 11.0  | 23.2   | 57.2   |  | StdDev:                        | 7.9   | 26.8   | 60.5   |      |
| % CV:                          | 6.2   | 4.1    | 4.5    |  | % CV:                          | 4.0   | 5.0    | 5.3    |      |

### Accuracy assessment

Comparison of LC-MS/MS quantitation of DHT in 54 donor serum samples between our research lab and to a reference lab is summarized in Table 2. DHT values ranged from 47 to 973 pg/mL. Only 4 out of 54 results differed by more than 20% and none were more than 22.4%. The differences between the two methods averaged 2.5%, a small positive bias by our lab. The two DHT methods were equivalent within an allowable total error (TEa) of 25%.

Table 2. Comparison of donor serum DHT results (pg/mL) between reference lab (Ref Lab) and research lab (BRL).

| Sample | Ref Lab | BRL | % Diff | Sample | Ref Lab | BRL  | % Diff | Sample | Ref Lab | BRL | % Diff |
|--------|---------|-----|--------|--------|---------|------|--------|--------|---------|-----|--------|
| 1      | 151     | 151 | 0.00   | 19     | 182     | 182  | 0.00   | 37     | 937     | 891 | -4.91  |
| 2      | 160     | 191 | 19.38  | 20     | 535     | 521  | -2.62  | 38     | 174     | 213 | 22.41  |
| 3      | 72      | 76  | 5.56   | 21     | 1380    | 1492 | 8.12   | 39     | 552     | 560 | 1.45   |
| 4      | 91      | 105 | 15.38  | 22     | 286     | 330  | 15.38  | 40     | 285     | 284 | -0.35  |
| 5      | 685     | 663 | -3.21  | 23     | 380     | 452  | 18.95  | 41     | 34      | 35  | 2.94   |
| 6      | 281     | 264 | -6.05  | 24     | 777     | 752  | -3.22  | 42     | 122     | 142 | 16.39  |
| 7      | 92      | 105 | 14.13  | 25     | 456     | 489  | 7.24   | 43     | 139     | 149 | 7.19   |
| 8      | 290     | 295 | 1.72   | 26     | 266     | 232  | -12.78 | 44     | 211     | 227 | 7.58   |
| 9      | 572     | 537 | -6.12  | 27     | 1070    | 1080 | 0.93   | 45     | 319     | 278 | -12.85 |
| 10     | 47      | 55  | 17.02  | 28     | 252     | 255  | 1.19   | 46     | 457     | 504 | 10.28  |
| 11     | 157     | 148 | -5.73  | 29     | 691     | 589  | -14.76 | 47     | 195     | 171 | -12.31 |
| 12     | 77      | 80  | 3.90   | 30     | 1160    | 1240 | 6.90   | 48     | 117     | 136 | 16.24  |
| 13     | 107     | 130 | 21.50  | 31     | 1230    | 1160 | -5.69  | 49     | 279     | 333 | 19.35  |
| 14     | 326     | 376 | 15.34  | 32     | 1020    | 804  | -21.18 | 50     | 259     | 239 | -7.72  |
| 15     | 462     | 408 | -11.69 | 33     | 55      | 44   | -20.00 | 51     | 470     | 480 | 2.13   |
| 16     | 76      | 69  | -9.21  | 34     | 255     | 240  | -5.88  | 52     | 449     | 493 | 9.80   |
| 17     | 66      | 75  | 13.64  | 35     | 339     | 311  | -8.26  | 53     | 240     | 261 | 8.75   |
| 18     | 293     | 305 | 4.10   | 36     | 94      | 75   | -20.21 | 54     | 267     | 306 | 14.61  |

### Throughput

Single-channel throughput was 12 injections per hour. When multi-channelled across 2, 3 or 4 channels, the throughput increased to 24, 36 and 48 injections per hour, respectively. DHT batches were also multi-channelled with pregnenolone and estrogen batches, which utilized the same HESI source conditions.

## CONCLUSIONS

Dihydrotestosterone (DHT) can be accurately measured in blood serum by this research method to achieve:

- Analytical range from 25 to 2,000 pg/mL
- Inter- & intra-batch precisions less than 7% CV and carryover less than 0.2%
- Results virtually equivalent to reference lab results
- Throughputs of 12, 24 or 48 injections per hour from a 1-, 2- or 4-channel system
- Multi-channeling with other methods utilizing the same HESI source conditions

## REFERENCES

- J. DiBussolo & M. Kozak, High-Throughput LC-MS/MS Quantification of Pregnenolone in Human Blood Serum for Clinical Research Purposes, Thermo Scientific Application Note AN-650, 2016.
- J. DiBussolo, How to Electro-Clean APCI Sources Between Injections, Thermo Scientific Poster Note PN64408, 2016.

## TRADEMARKS/LICENSING

© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

MSACL EU 2018 Poster

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.